

# **Instruction and Operation Manual**

**S605** 

# **Portable Breathing Air Quality Analyzer**



## .SUO

Dear Customer,

Thank you for choosing our product.

Before starting up the device, please read the manual in full and observe the operating instructions stated. The manufacturer cannot be held liable for any damage which occurs as a result of non-observance or noncompliance with this manual.

Should the device be tampered with in any manner other than a procedure which is described and specified in the manual, the warranty is void and the manufacturer is exempt from liability.

The device is designed exclusively for the described application.

We offer no guarantee for the suitability for any other purpose. We are also not liable for consequential damage resulting from the delivery, capability or use of this device.

2 \$605

# **Table of Contents**

1	Safety Instructions	5
2	Registered Trademarks	8
3	Application	8
4	Features	8
5	Technical Data	9
	5.1 General Data	9
	5.2 Operating Conditions	
	5.3 Measurement Data	
6	Dimensional Drawings	
	Installation on Site	
	7.1 Compressed Air Inlet	
	7.2 Compressed Air Connections	
	7.2.1 Installation Requirements	
	7.2.2 Installation	
	7.3 Electrical Connectors	.17
	7.4 Electrical Connections	.17
	7.4.1 RS-485 Networking (Modbus/RTU)	
	7.4.2 TCP/IP Networking (Modbus/TCP)	
	7.4.3 Connect with a PC or an OTG USB Stick	
	7.5 Water Filling Connector	
	7.6 Calibration Connector	
8	Setup and Configuration	.20
	Operation	
	9.1 Main Screen	
	9.1.1 Quick Buttons and Icons	
	9.1.2 Icons in the Status Bar	
	9.2 Value View	.24
	9.3 Graphic View	
	9.4 Menu	
	9.5 Configure S4A Remote Connection	
	9.6 SIM Card Requirement for 4G/LTE Option	
	9.6.1 Data Traffic	
	9.6.2 PIN Code Protection	
1	0 Guided Measurement	
_	10.1 Steps for Guided Measurement	
	10.2 Reports for Guided Measurements	
	•	

## .SU ()

11 Firmware Updates	33
12 Calibration	35
12.1 Calibration for DP Sensor and Oil Sensor	35
12.2 Calibration for O2, CO, CO2 Sensor	35
12.2.1 Calibration Gases	
12.2.2 Calibration Procedure	36
13 Optional Accessories	42
14 Maintenance	42
15 Disposal or Waste	42
16 Warranty	42
17 Appendix - Modbus Holding Register Table	
17.1 Modbus Interface	43
17.2 Modbus Register Table	44

4 \$605

## 1 Safety Instructions



## Please check if this instruction manual matches the product type.

Please observe all notes and instructions indicated in this manual. It contains essential information which must be observed before and during installation, operation and maintenance. Therefore this instruction manual must be read carefully by the technician as well as by the responsible user / qualified personnel.

This instruction manual must be available at the operation site of the flow sensor at any time. In case of any obscurities or questions, regarding this manual or the product, please contact the manufacturer.



#### **WARNING!**

## Compressed air!

Any contact with quickly escaping air or bursting parts of the compressed air system can lead to serious injuries or even death!

- Do not exceed the maximum permitted pressure range (see sensors label and/or the specifications written in this manual).
- Only use pressure tight installation material.
- Avoid that persons get hit by escaping air or bursting parts of the instrument.
- The system must be pressure less during maintenance work.



#### **WARNING!**

## Voltage used for supply!

Any contact with energized parts of the product may lead to an electrical shock which can lead to serious injuries or even death!

- Consider all regulations for electrical installations.
- The system must be disconnected from any power supply during maintenance work.
- Do not exceed the permitted operating parameters as specified in the data sheet, this manual or written on the instrument.
- Make sure the product is operated in its permitted limitations.



- Do not exceed or undercut the permitted storage and operation temperature, pressure and humidity.
- The product should be maintained and calibrated frequently, at least annually. For more information, please contact the customer support.

## **General safety instructions**

- It is not allowed to use the product in explosive areas.
- Please observe the national regulations before/during installation and operation.

#### Remarks

- It is not allowed to disassemble the product.
- Always check the compressed air connectors in terms of stability and tightness before any operation.



#### **ATTENTION!**

Measurement values can be affected by malfunction!

The product must be installed properly and frequently maintained, otherwise it may lead to wrong measurements, which can lead to wrong results.

- The instrument must be operated within or better than compressed air quality Class 4.4.4 according to the ISO8573-1. If the air quality cannot reach the specified class, the device may show wrong measurement values or even may get damaged.
- Do not exceed the maximum operation temperature or pressure.
- Avoid condensation inside the instrument caused by the supplied air or gas because it harms the instrument and affects the accuracy.
- It is recommended that before any operation you examine the instrument in terms of the measurement accuracy and deviations. To do this, you may use high efficiency filters for oil vapor and particle, and the instrument should show low values close to zero. Additionally, it is recommended you examine the dew point measurement to ensure it is working correctly. To do this, you may use a dry gas with a defined dew point value which should be reached within several minutes. Or you may use any reference meter. Only after doing these examinations, you can rely on the measurement results after the instrument have been used for other measurement cycles.



 Always perform a visual inspection in prior to every measurement before connecting compressed air. If any loose objects are found, do not connect the compressed air because it may lead to a damage or injuries.

### Storage and transportation

- Make sure that the transportation temperature is between -10... 50°C.
- For transportation it is recommended to use the packaging which comes with the sensor.
- Please make sure that the storage temperature of the sensor is between -10... 50°C.
- Avoid direct UV and solar radiation during storage.
- For the storage the humidity must be < 90%, no condensation.




#### **ATTENTION!**

## **Equipment may get damaged!**

Please make sure, that your measuring point is free of excessive contamination/dirt. This should maintained before every measurement.

- Observe the measuring point always before measurement if it is free of contamination like water drops, oil drops or other rough contaminants.
- Should water hit the inner electronics, the senors could be seriously damaged.
- Any rough contamination of particles, water, oil or other impurities may damage the instrument or affect the accuracy.
- Check your measurement point with the enclosed test kit, by purging air in prior to the use of the instrument through the test kit. A clean tissue may help to identify oil, water or any other lubricants in the air.





## 2 Registered Trademarks

SUTO<sup>®</sup> Registered trademark of SUTO iTEC

MODBUS<sup>®</sup> Registered trademark of the Modbus Organization,

Hopkinton, USA

## 3 Application

The S605 Portable Breathing Air Quality Analyzer measures the content of  $O_2$ ,  $CO_2$ ,  $H_2O$  and Oil vapor according to the breathing air purity standards. It monitors the contaminant concentration to ensure that the breathing air quality complies with international standards such as EN 12021 or CFSR 1910.134(d).

The S605 is mainly used in an industrial environment. It can not be used in hazardous or explosive environment and it not suited for permanent outside applications.

#### 4 Features

- All-in-one device measures O<sub>2</sub>, CO, CO<sub>2</sub>, H<sub>2</sub>O, Oil, and pressure.
- Portable design makes it to be carried with one hand.
- The device has a DC power supply and additionally a battery which can bridge power failures.
- High resolution 5" color touchscreen display.
- Software guided measurement.
- USB port for transferring data to a USB-OTG memory.
- Modbus/TCP and Modbus/RTU interfaces for communication.
- Data logger with the capacity of 30-million values.
- Only one gas inlet for all parameters.
- 0.4...1.5 MPa(g) pressure at inlet; higher pressure is available with external pressure reducer.
- With optional 4G modem feature (P/N: A1670) and dedicated S4A software, you can check the online measurement values, read out the logging files remotely.



## **5 Technical Data**

## 5.1 General Data

CE FC CA		
Display	Touchscreen, Size: 5", Resolution: 800 x 480 pix	
Power supply	Mains supply adapter (AC/DC) Input: 100 240 VAC, 50/60 Hz, 1.8 A Output: 20 VDC, 3.25 A, 65 W max	
Internal battery	7500 mAh Rechargeable for nearly 3.5 hours operation Charging time about 1.5 hours	
Process connection	Micro quick connector, full pass through, male (1.5 m hose with coupling included)	
Interface	USB M12 (Modbus/RTU) Ethernet (Modbus/TCP) USB-A for 4G/LTE modem (optional)	
Data logger	Internal, 30-million values	
Housing material	PC + ABS, Aluminum	
Protection class	IP54 (cover closed)	
Dimension	470 x 365 x 181 mm	
Approvals	FCC, EC	
Weight	11 kg	

## **5.2 Operating Conditions**

Measuring medium	Compressed breathing air
Humidity of the medium	< 40%, non condensation
Temperature of the medium	0°C +45°C
Operation pressure	0.4 1.5 MPa(g) Higher pressure needs an external pressure reducer.
Storage & Transport temperature	-10°C + 50°C
Ambient temperature	0°C +50°C



Ambient humidity	0 90% rH
Sampling flow rate	6 LPM@4 MPa(g), depends on input pressure
Sampling rate	1 sample/sec

## 5.3 Measurement Data

## Oil Vapor Sensor (only for S605-I version)

Measuring range	0.001 5.000 mg/m³ (Based on 1000 hPa(a), 20°C, 0% relative humidity)
Accuracy	5% of reading $\pm$ 0.003 mg/m <sup>3</sup>
Resolution	0.001 mg/m <sup>3</sup>
Sensor principle	Photo ionization detector
UV lamp lifetime	6,000 working hours or 1 year, whichever comes first

## Oil Mist and Particle Sensor (only for S605-C version)

Measuring range	0.0 5.0 mg/m³ (Based on 1000 hPa(a), 20°C, 0% relative humidity)
Accuracy	15% of reading ± 0.1 mg/m <sup>3</sup>
Resolution	0.1 mg/m <sup>3</sup>
Sensor principle	Oil mist and particle sensor

## O<sub>2</sub> Sensor

Measuring range	0 25%	
Accuracy	±(1% of reading +0.05%)	
Resolution	0.1%	
Sensor principle	Optical oxygen sensor	
CO Sensor		
Measuring range 0 20 ppm		
Accuracy	±(5% of reading +1 ppm)	
Resolution	0.1 ppm	
Sensor principle Electrochemical sensor		
CO₂ Sensor		

10 \$605



Measuring range	0 1000 ppm		
	<u> </u>		
Accuracy	±(1% of reading +25 ppm)		
Resolution	1 ppm		
Sensor principle	NDIR sensor		
H <sub>2</sub> O			
Measuring range	-100 +20 °C Td / 0 17458.6 mg/m <sup>3</sup>		
Accuracy	±1 °C Td (0 20 °C Td) ±2 °C Td (-700 °C Td) ±3 °C Td (-10070 °C Td)		
Resolution	0.1 °C Td		
Sensor principle	QCM + Polymer		
Pressure Sensor			
Measuring range 0 16 bar(g)			
Accuracy	0.5% FS		
Resolution	0.01 bar		
Sensor principle	Piezzo resistive pressure sensor		



## **6 Dimensional Drawings**

Dimension of S605 in mm (cover closed):



12 \$605



## 7 Installation on Site

Please make sure that all components listed below are included in your package.


Qty.	Description	Item No.	
1	S605 Portable Breathing Air Quality Analyzer, in a hand carrying case with handle and shoulder belt		
1	USB OTG memory stick	S605-I version:	
1	Purge filter for pre-measurement (test kit)	P560 0605 (with oil vapor sensor)	
1	Power adapter Input: 100 240 VAC,50/60 Hz, 1.8A Output: 20 VDC, 3.25 A, 65 W max	S605-C version: P560 1605 (with	
1	1.5 m Teflon hose with a quick connector and a compressed air coupling at each end	oil mist and particle sensor)	
1	Water filling bottle		
1	Operation and instruction manual		
1	Certificate of calibration	-	
1	M12 connector		
1	USB 4G/LTE modem, including S4A software (option)	A1670	

To replace the items in the preceding list or buy accessories, contact the manufacturer or your local distributor.



## 7.1 Compressed Air Inlet

The compressed air inlet is located on the right side of the S605 housing.



S605-I right view



#### **ATTENTION!**

**Permissible pressure!** 

Please observe the maximum permissible inlet pressure. If the pressure exceeds the range, the device will be damaged. If the pressure is too low, the volumetric flow will not be high enough, which leads to wrong results.



## 7.2 Compressed Air Connections

### 7.2.1 Installation Requirements

Make sure that the installation meets the following requirements:

- The device must be set up next to the measuring point. Make sure that the device is put on a flat place.
- The tubes should be not bended heavily and should be installed with a big curve radius to avoid turbulence in the air flow.
- Connect the S605 to the power during the measurement, and make sure that the device is not turned off or plugged off during the measurement because data will then be lost.



#### **ATTENTION!**


If the device is not installed properly it may lead to wrong measurement results.

- The device is designed to be operated indoors only. To use the device outdoors, protect the device from direct sunlight and rain.
- Before you start the measurement, check your point of measurement for contaminants.



#### **ATTENTION!**

Before you connect the device to your point of measurement, make sure that there is no rough contamination like water/oil drops or heavy dust. This may damage the sensor units. For this please use the purge filter test kit.





Before each measurement, the measuring point should be checked for coarse contamination. Follow the steps below:

- 1. Connect the purge filter test kit onto your measuring point first. Turn on the purge valve on the test kit and purge air for a short period.
- 2. Check the filter in the test kit to identify whether it shows high contamination of water, oil or dust.
- 3. If the filter is contaminated severely, stop using the S605 for measurement because this may lead to serious damage to the device. In case you are not sure, please contact the manufacturer.

#### 7.2.2 Installation

- 1. Remove the protection cap from the air inlet of the S605.
- 2. Connect your compressed air system to the air inlet of the S605 by using the supplied Teflon hose.


#### **Notes:**

- Make sure you are using clean and oil-free components when connecting the instrument.
- The compressed air supply must always be applied carefully. If you
  are using ball valves, turn on them carefully and slowly. If you
  observe any unexpected noise or air leakage, shut down the
  compressed air supply immediately and disconnect the instrument.



#### 7.3 Electrical Connectors

Three connectors are available on the right side of the S605.



S605-I right view

Туре	Description
RJ-45	Modbus/TCP
M12	Modbus/RTU
Type-C USB	Power supply connector /Battery charger interface



#### **ATTENTION!**


Please connect the power correctly! If wrong supply is used such as AC power, the instrument may get damaged.

## 7.4 Electrical Connections

## 7.4.1 RS-485 Networking (Modbus/RTU)

Through the M12 connector, the S605 can be connected to the RS-485 network over the Modbus/RTU protocol.

The RS-485 networking uses the following 3 pins in the M12 connector. The pin functions are described below.



Pin No.	Assignment	Function
4	Data +	Data + signal
5	Data -	Data - signal
1	GND RS-485	Ground of RS-485



## 7.4.2 TCP/IP Networking (Modbus/TCP)

Through the RJ-45 connector, the S605 can be connected to the TCP/IP network over the Modbus/TCP protocol.

Remove the protection cap and plug in the network cable (RJ-45).

For more information about the data attributes of measurement channels, see Appendix - Modbus Holding Register Table.

#### 7.4.3 Connect with a PC or an OTG USB Stick


Through the USB port on the front panel, the S605 can be connected with:

- An OTG memory stick: To import firmware for upgrade and to export data.
- A PC where data analysis software such as S4A or S4M is installed: To export data to the PC for analysis.

## 7.5 Water Filling Connector

To ensure that the measured gas is in a reasonable humidity, a humidifier is integrated in the S605, thus the CO sensor will not be damaged because the measured gas is too dry.

The water filling connector of the humidifier is on the front side of the S605. The water draining connector (i.e., air outlet) is on the right of the S605.



S605-I right view

When the water level is too low, the low water level indication in the status bar blinks, which reminds to fill water for the humidifier.

A water filling bottle is supplied in the package. To fill water, remove the

protection caps from the water inlet and air outlet first, plug bottle nozzle into the water filling connector, then squeeze the bottle and the water is filled in. This process may take 1-2 minutes. Obersve the water level indication in the status bar. When the high water level icon appears on the status bar, stop to fill. Then put on the two protection caps.



#### **ATTENTION!**

Please fill with water before the first use!


After the S605 is functional, refill water in time.

The humidifier lifetime is nearly 18 days (base on the general condition). If the water level is too low, it will greatly reduce the lifetime of the CO sensor or even damage the sensor!

For the injected water, it is recommended to use pure/demineralized water. Unclean water can cause scaling and even contaminate sensors.

#### 7.6 Calibration Connector

The calibration connector is located on the right side of the S605. When the sensors need to be calibrated, remove the connector cap, then connect the calibration hose. When the calibration is finished, put on the cap to avoid the dust or other substance entering. For more information on calibration procedure, see section 12 Calibration.



S605-I right view



## 8 Setup and Configuration

The S605 is configured ex-works and ready to work when you get it. The S605 provides a guided measurement procedure to take you through device setup for each measurement parameter. All these setup settings are automatically saved into the device even on a power failure.

For more information about the guided measurement procedure, see section 10 Guided Measurement.



#### **ATTENTION**

If you are facing problems setting up your device, contact the manufacturer or your local dealer for assistance.

If power outage occurs during the measurement, the measurement data will not be saved!




## 9 Operation




After the S605 is powered on, the initialization screen is displayed with an active progress bar.

During the initialization, the device configures the sensors and runs initialization routines.




After the initialization is completed, the system information is shown including device type, serial number, and more. This information displays for 3 seconds then moves to the next page.




The sensor information screen shows the important maintenance information of the sensor. This screen displays for 3 seconds then moves to the next page.





The attention screen reminds to use the test kit to avoid the sensor damage. Click **OK** to move to the next screen.




This screen asks if need to input the system pressure.

If you perform a high

If you perform a high pressure measurement (>15 barg) through a pressure reducer, enter the system pressure before the pressure reducer.

Click **Yes** to input the system pressure. Clicking **No** goes to the Value screen.

**Note**: If your system pressure is greater than 15 bar, be sure to enter your system pressure correctly, otherwise it may affect the measurement!




The value screen appears, as shown on the left.



#### 9.1 Main Screen

The following figure shows the main screen in the value view.



Status bar

Display area

Quick buttons and icons

## 9.1.1 Quick Buttons and Icons

<b>Buttons and icons</b>	Description
Graphic	To switch to the graphic view of measurement data
Value	To switch to the value view of measurement data
Menu	To access the operation menus
Ō	To take a screenshot

#### 9.1.2 Icons in the Status Bar



System errors occur. Press this icon to get further information.



Calibration is expired.

Contact the manufacturer or your local dealer.



Low water level indication



High water level indication



4G/LTE signal strength



S4A remote connected



Data logger status:


- STOP—Indicates that the data logger is not running.
- LOG—Indicates that the data logger is running.



#### 9.2 Value View

All measured values, channel range, limit value, alarms can be viewed in realtime.

The green color bar area indicates the valid range; the red color bar area indicates the invalid range. When the measured value exceeds the alarm threshold, the data color will be red and keep blinking.



To switch to the value view, click **Value** in the bottom.

## 9.3 Graphic View

The measured result is shown in a curve when the recording is activated. It shows the last 30 minutes by default.



To switch to the Graphic view, click **Graphic** on the bottom.

The graphic view is pre-configured in the factory, and you do not need to change anything.


If you make changes, follow the description below.

Item	Description
Y-axes	Touch Y-axes to scale it.
X-axes	Touch X-axes to define viewing period.



Item	Description
< >	Scroll time.
<b>♠</b>	Click the icon to go back to the current time.
Channel data area	<ul><li>Short touch can switch through single channel.</li><li>Long touch is for settings.</li></ul>

#### 9.4 Menu



To access the main menu, click **Menu** in the bottom bar.

All operation menus are displayed.

These operation menus are described below.

#### **Guided Measurement**

To start the guided measurements, which lead you through a complete measurement cycle.

For more information, see Chapter 10 Guided Measurement.

## **Sensor settings**

To choose the different standards. When choosing the standard, the limit value is automatically set. The customer also can manually set the limit value according to requirements.

## System settings

To configure system-level settings:

- To perform general settings on access code, date, time and language and so on.
- To view device information, such as the serial number, firmware and hardware versions.
- To calibrate the touch screen if it does not respond to user inputs correctly or precisely.
- To update the system firmware.



#### **Service Company info**

Show the service company contact information. You can import a logo file to the S605 through the USB port.

#### **Files**

To copy and delete the report files, recording files, calibration record files, captured screenshots and check the memory status.

#### **Communication**

To configure communication parameters. Only the following sub-menus are applicable to S605:


- Modbus/TCP: The S605 can be connected to the TCP/IP network over the Modbus/TCP protocol. In this sub-menu, select the communication protocol as required, and configure the IP address of the S605.
- **Modbus/RTU**: The S605 can be connected to the RS-485 network over the Modbus/RTU protocol. You can configure the communication parameters such as address, baud rate etc.
- **S4A Remote**: Used to configure the S605 to connect to the S4A remote server, viewing the measurement values, and reading out the logging files via the S4A software.

#### **Calibration**

This is used for the zero and span calibration for the  $O_2$ , CO,  $CO_2$  sensors. You can check the calibration information for all the sensors include the exchange sensor unit.



## 9.5 Configure S4A Remote Connection



Follow steps below to connect S4A remote:

1. Install a SIM card into the provided 4G/LTE Modem and insert the 4G/LTE modem into the USB port.

**Note**: Turn off the S605 before inserting the 4G module, and turn on the S605 after the 4G module is inserted.

- 2. Click **Menu** > **Communication** > **S4A** remote to enter the S4A configuration page.
- 3. Input the device ID and password. Click **Connect** to connect the S605 to the S4A remote server.
- 4. Start the S4A software to establish the connection. Then you can view the online measurement values, and read out the logging files remotely with the S4A software.

For the detailed configuration steps, watch the instruction video which can be found on SUTO website.

## 9.6 SIM Card Requirement for 4G/LTE Option

To establish the 4G/LTE communication between S605 and S4A software with S4A remote function, customer needs to buy a SIM card locally.

#### 9.6.1 Data Traffic

Data traffic required is 7 GB/month if the system is running  $7 \times 24$  hours.

#### 9.6.2 PIN Code Protection

The S605 does not support modifying PIN code on the screen. If your




SIM card has PIN code protection, disable it by using other devices, for instance using your mobile phone.

#### 10 Guided Measurement


The S605 provides a software-based guided measurement which takes you through the complete measurement. This leads to a simplified measurement process and prevents you from wrong measurements.

To start a guided measurement, do the following:

1. Click Menu > Guided measurement.



2. Select the type of measurement that you want to perform.






- Monitoring with programmed time out: It starts a
  measurement with a user-programmed period of measurement
  time. You can set the measurement time during the process of
  measurement preparation. The system will then, after finishing the
  programmed measurement duration, stop the measurement
  automatically and save the data. This mode is ideally used for
  audits where you must measure at several points. You can program
  for each point a duration of e.g. 2 hours and then you can compare
  the measurements.
- **Monitoring with manual stop**: It starts the measurement without a programmed stop time. You can click it to start the measurement and stop it whenever you want. This mode can be used to monitor changes of values.
- **Monitoring until values valid:** The measurement will stop when all values are within valid ranges for a certain time. This mode records the complete alarm information. The user can freely set the time out according to the requirement.
- 3. To perform the guided measurement, follow the onscreen instructions. For more information, see section 10.1 Steps for Guided Measurement.
- 4. Before starting the measurement, it is recommended to select the standard to be followed first. The S605 will set the alarm threshold according to the selected standard, evaluate the air quality according to the standard, and show all the information in the report.
- To manage the measurement files generated, click Menu > Files > Reports.



## 10.1 Steps for Guided Measurement

After you start a guided measurement, follow steps below to go through the whole process.




 Input the location information (that is, customer's data), which will be shown in the report.

Enter how long the measurement will last.
 The longer it lasts, the more stabilized the values will be and the more closer it will approximate to the system reality.

#### Notes:

- In order to ensure the best measurement results, it is recommended to set the measurement time greater than 45 minutes.
- The step 2 comes up only when the Monitoring with programmed time out is selected.
- The Monitoring with manual stop indicates that no measurement duration needs to be entered.
- The **Monitoring until values valid** indicates the measurement duration and the valid time duration need to be entered.






3. A confirmation displays before the measurement starts. Click **Yes** to start.

4. The system checks whether the pressure is in the valid range or not.

5. A pop-up appears and asks if the cleaning is needed. If you click **Yes**, the cleaning process starts. If you click **NO**, the measurement starts.






During the measurement, the data logger status icon on the status bar switches from **STOP** to **LOG**.

The remaining time is displayed on the upper left corner.


Please wait until the system completes the measurement.



When the measurement is successfully completed, the screen on the left appears. Click **OK** to finish the measurement.

## 10.2 Reports for Guided Measurements

After performing guided measurements, you can view and manage measurement files through **Menu** > **Files** > **Reports**.



To copy or delete files, select the file radio boxes, and click the corresponding button at the bottom.




## 11 Firmware Updates

To update the firmware, just need a OTG-USB drive(supplied with the S605) and the firmware file (format like DIS605\_1.82.tar)




1. Insert the USB drive to the connector on the front panel.

Make sure to do this after the S605 has fully booted.




2. Click **Menu** to enter the Menu page.




3. Click **System** settings.





4. Click **System Update**.



5. A pop up will be shown. Then click on the firmware file in this example DIS605\_1.82.tar.



6. A new pop up will ask you to confirm the update. Click **Yes** to start the update. Click **NO** to abort the operation.

#### Notes:

Make sure you do not turn off the power.

Make sure you do not disconnect the USB drive.


Wait until the device is fully booted again before removing the USB drive.

After the update is done, you can select the new language or use new functions offered in the update.



#### 12 Calibration

Because the sensor accuracy is affected by the environment conditions, it is recommend to calibrate the device at least once per year. The calibration should be done by specialized engineers. You can check calibration information through **Menu** > **Calibration** > **Calibration Information**.



#### 12.1 Calibration for DP Sensor and Oil Sensor

For the dew point sensor and oil sensor, the calibration operation is complicated and time consuming, so SUTO provides the exchange service, that is, the complete sensor unit is exchanged with a newly-calibrated unit. This saves time and costs.

## 12.2 Calibration for O<sub>2</sub>, CO<sub>2</sub> CO<sub>2</sub> Sensor

The  $O_2$  sensor, CO sensor, and  $CO_2$  sensor can be calibrated on site or replaced with calibrated sensors using the exchange service of SUTO.



#### **ATTENTION!**

Please save all your measurement data on an external device before returning the instrument to calibration and service. It might be necessary to reset the displays storage during calibration and service.



#### 12.2.1 Calibration Gases

Depending on the type of calibration, different calibration gases are used.

#### Zero gas

The zero gas is a test gas used for zero-point adjustment. The zero gas must be free from the measured gas or other substances that have cross-sensitivity.

### Calibration and test gas

The span gas is a calibration gas used for sensitivity adjustment or verification of the calibration result. The span gas used can be the measured gas.

Here is recommended gases for the calibration and test gases:

Zero gas: 8%...12% O<sub>2</sub>, N<sub>2</sub> background

Calibration gas O<sub>2</sub>: 20%...24% O<sub>2</sub>, N<sub>2</sub> background

Calibration gas CO: 16...20 ppm CO, N<sub>2</sub> background

Calibration gas CO<sub>2</sub>: 800...1000 ppm CO<sub>2</sub>, N<sub>2</sub> background

The test gas can be the ones in the middle. For example:

Test gas CO: 10 ppm CO,  $N_2$  background(the calibration gas is 20 ppm CO)

Test gas  $CO_2$ : 500 ppm  $CO_2$ ,  $N_2$  background(the calibration gas is 1000 ppm  $CO_2$ )

Uncertainty of gases is 1% (K=2)




#### Attention!

Please select the suitable calibration gas, otherwise the sensor measurement accuracy will be affected.

#### 12.2.2 Calibration Procedure

You can enter the calibration through **Menu** > **Calibration**.








#### **WARNING!**

Do the calibration after the sensor is ready to work, for example the warm-up phase of the sensor is completed.

Do not inhale the calibration gas. It can be harmful to health or even death!

#### **User Calibration Settings**

The calibration reference data can be set before the calibration through **Menu** > **Calibration** > **User Calibration Settings**.




#### **Zero calibration**

Enter the zero calibration through **Menu** > **Calibration** > **Zero** 




#### calibration.



#### Span calibration

Enter the span calibration through **Menu** > **Calibration** > **Span calibration**.



### Zero/Span calibration procedures

- 1. Cut off air supply from main inlet.
- 2. Screw off the cap from calibration connector if not.
- 3. Connect the calibration hose from the gas cylinder to the calibration connector.
- 4. Adjust the flow controller at the gas cylinder to 0.5 l/min and the output pressure must 0 bar(g).
- 5. Select Menu > Calibration > Zero/Span calibration.




- 6. Select the sensor to be calibrated.
- 7. Supply zero/span gas. Make sure the flow rate of 0.5L/min.
- 8. Click **Start** button, and a progress bar will show with the calibration remaining time.
- 9. Zero/Span calibration takes nearly ten minutes. When the calibration is finished, you can check the calibration result in the **Status** field.
- 10. After the calibration is done, close the gas supply and remove the hose.

After the calibration, you can check the calibration result, by using another span gas. Just perform the same steps as above, but select the **test** in the span calibration menu.


#### **Calibration Records**

You can check the calibration records through **Menu** > **Files** > **Calibration Records**.




Example for calibration records file:





The figure below is an example for the gas connection.





## Here is the recommended supplier and accessories list:

Item	Brand	Description	Picture
Flow controller	MESA	Series 400 Fixed Flow Regulator ID:400-LPM	80 1 1 mm Park M 10 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1
Calibration Gas	MESA	Zero gas O <sub>2</sub> : 11% O <sub>2</sub> , N <sub>2</sub> background O <sub>2</sub> : 22% O <sub>2</sub> , N <sub>2</sub> background CO: 20 ppm CO, N <sub>2</sub> background CO <sub>2</sub> : 1000 ppm CO <sub>2</sub> , N <sub>2</sub> background 34L, steel calibration gas cylinders	CHERRATION AND SECURITY AND COMPANIES AND CO
Test Gas	MESA	CO: 10 ppm CO, N <sub>2</sub> background CO <sub>2</sub> : 500 ppm CO <sub>2</sub> , N <sub>2</sub> background 34L, steel calibration gas cylinders	CLEGANO DE SIA LUCIO DE COMPOSITO DE COMPOSI
Transport casing	MESA	Carrying cases ID:500	

Website: <a href="https://mesagas.com/">https://mesagas.com/</a>



### 13 Optional Accessories

To purchase optional accessories, please contact the manufacturer or your local dealer:

- Teflon hoses and sorts of adapters
- USB 4G modem, including S4A software (P/N: A1670)
- Pressure reducer, including transport case (P/N: A604 0004)

#### 14 Maintenance

Use a moist fabric to clean the device. For the use in GMP areas, the device must be disinfected through wipe disinfection. For more information, please contact the manufacturer or your local dealer.



#### **ATTENTION!**

Please dry the device after cleaning using a clean and dry fabric. Always take care, that the fabric for cleaning is not wet because water could get into the device and lead to damage.

### 15 Disposal or Waste



Electronic devices are recyclable material and do not belong in the household waste.

The sensor, the accessories and its packing must be disposed according to your local statutory requirements. The dispose can also be carried by the manufacturer of the product, for this please contact the manufacturer.

### 16 Warranty

Please find the warranty as a separated warranty card included with the instrument delivery.

The warranty does not cover any wear parts or consumables, therefore the UV lamp with limited lifetime as well as the internal filter are not covered by the warranty.



# 17 Appendix - Modbus Holding Register Table

#### 17.1 Modbus Interface

The default settings of the Modbus interface are as follows:

#### Communication parameters (Modbus/RTU)

Baud rate : 19200

Device address : Last two digits of serial number

Framing / parity / stop bit: 8, N, 1

Response time : 1 second

Response delay : 0 ms
Inter-frame spacing : 7 char

#### **Communication parameters (Modbus/TCP)**

DHCP : Yes

MAC : Set ex-factory

IP address : Dynamic or Static Subnet : Dynamic or Static Gateway : Dynamic or Static

Timeout :  $\geq$  200 ms

### Response message that the device returns to the master:

Function code: 03

The information of the byte order is shown in the table below:

Byte	Sequence						
Order	1st	2nd	3rd	4th	Туре		
1-0-3-2	Byte 1 (MMMMMMM*)	Byte 0 (MMMMMMM *)	Byte 3 (SEEEEEEE)	Byte 2 (EMMMMMMM *)	FLOAT		
1-0-3-2	Byte 1	Byte 0 LSB	Byte 3 MSB	Byte 2	UINT32 INT32		
1-0	Byte 1 MSB	Byte 0 LSB			UINT16 INT16		
1-0	Byte 1 XXX *	Byte 0 DATA			UINT8 INT8		

<sup>\*</sup> S: Sign, E: Exponent, M: Mantissa, XXX: no value



#### **Explanations of MSB and LSB**

MSB MSB refers to Most Significant Byte first, which follows the Big-Endian byte order.

For example, if the main system follows the MSB first order: When the 4-byte floating number, in the order of Byte1-Byte0-Byte3-Byte2, is received from the slave (sensor), the master must change the byte order to Byte3-Byte2-Byte1-Byte0 for the correct display of the value.

LSB refers to Least Significant Byte first, which follows the Little-Endian byte order.

For example, if the main system follows the LSB first order: When the 4-byte floating number, in the order of Byte1-Byte0-Byte3-Byte2, is received from the slave (sensor), the master must change the byte order to Byte0-Byte1-Byte2-Byte3 for the correct display of the value.

#### 17.2 Modbus Register Table

The following table lists specifications of the Modbus output channels in this device.

Holding register	Data type	Data Length	Channel Description	R/ W	Remark			
	Common System Information							
2000	INT16U	2-Byte	Group ID (=6)	R	Gas flow sensor, water flow etc.			
2001	INT16U	2-Byte	Device ID	R	Unique identification			
2002	INT32U	4-Byte	Serial number	R	Factory S/N			
2004	INT16U	2-Byte	High = 0xFF, Low byte HW version	R	Before it was SW version in High Byte, but in future there will be a 3 BYTE version and the High Byte is set to 0xFF			
2005	DOUBLE	8-Byte	Calibration date	R	Date time (8-bytes) format is time- variant (Microsoft)			
2009	INT16U	2-Byte	Valid days from calibration date	R				
2010	INT16U	2-Byte	Number of Measuring Channels	R	196			
2011	STRING	16-Byte	Device description	R/W	"S605" etc.			



2019	INT16U	2-Byte	Version of Modbus Definition	R	0x0201 = 02.01
2020	INT16U	2-Byte	Version of this register table	R	Start from 0x100 =1.00
2021	INT32U	4-Byte	Calibration date	R	unix time format
2023	INT32U	4-Byte	Production date	R	unix time format
2025	INT32U	4-Byte	FW version (low 3 Bytes used)	R	New firmware version format: XX.XX.XX
2027	STRING	16-Byte	Device Location	R/W	User can enter a text
2035	STRING	16-Byte	Measurement Point	R/W	User can enter a text
2043	INT16U	2-Byte	Number of sensors	R	How many and which SUTO sensors are connected to the main unit. Their settings are found from 10000 up.
			Modbus/RTU Settings		
2100	INT16U	2-BYTE	Write Password	W	Response: $0 \rightarrow access$ permitted $1 \rightarrow authorization$ required
					password: 0XA5A5
2110	INT16U	2-Byte	Modbus baud rate index 01200 12400 24800 39600 419200 538400 657600 7115200	R	
2111	INT16U	2-Byte	Modbus parity	R	
2112	INT16U	2-Byte	Modbus stop bit	R	
2113	INT16U	2-Byte	Device Address	R	
			Modbus/TCP Settings		
2120	INT32U	4-Byte	Static IP address	R	
2122	INT16U	2-Byte	Portal number	R	
2123	INT32U	4-Byte	Subnet mask	R	
2125	INT16U	2-Byte	DHCP enabled	R	
2126	INT32U	4-Byte	IP address	R	For DHCP only
			Alarm Settings		



2130	FLOAT	4-Byte	O2 High alarm threshold	R	
2132	FLOAT	4-Byte	O2 Low alarm threshold	R	
2134	FLOAT	4-Byte	O2 Hysteresis	R	
					For S606: bit7-bit0 (low byte): 0x00= relay 0
					0x07 = relay 7
					relay 0-7 are on the optional relay module
2136	INT16U	2-Byte	Relay port/status	R	bit15: 0= channel alarm relay close 1= channel alarm relay open
					bit 14: 0 = summary relay close 1 = summary relay open
					For S605: Always 0XFFFF: don't support
2137	FLOAT	4-Byte	CO2 High alarm threshold	R	
2139	FLOAT	4-Byte	CO2 Low alarm threshold	R	
2141	FLOAT	4-Byte	CO2 Hysteresis	R	
2143	INT16U	2-Byte	Relay port/status	R	See explanation at 2136
2144	FLOAT	4-Byte	CO High alarm threshold	R	
2146	FLOAT	4-Byte	CO Low alarm threshold	R	
2148	FLOAT	4-Byte	CO Hysteresis	R	
2150	INT16U	2-Byte	Relay port/status	R	See explanation at 2136
2151	FLOAT	4-Byte	H2O High alarm threshold	R	This is the dew point measurement.
2153	FLOAT	4-Byte	H2O Low alarm threshold	R	
2155	FLOAT	4-Byte	H2O Hysteresis	R	
2157	INT16U	2-Byte	Relay port/status	R	See explanation at 2136
2158	FLOAT	4-Byte	Oil High alarm threshold	R	
2160	FLOAT	4-Byte	Oil Low alarm threshold	R	



2162	FLOAT	4-Byte	Oil Hysteresis	R	
2164	INT16U	2-Byte	Relay port/status	R	See explanation at 2136
2165	FLOAT	4-Byte	Pressure High alarm threshold	R	
2167	FLOAT	4-Byte	Pressure Low alarm threshold	R	
2169	FLOAT	4-Byte	Pressure Hysteresis	R	
2171	INT16U	2-Byte	Relay port/status	R	See explanation at 2136
2172	FLOAT	4-Byte	Temperature High alarm threshold	R	
2174	FLOAT	4-Byte	Temperature Low alarm threshold	R	
2176	FLOAT	4-Byte	Temperature Hysteresis	R	
2178	INT16U	2-Byte	Relay port/status	R	See explanation at 2136
2179	FLOAT	4-Byte	Humidity to gas sensor High alarm threshold	R	
2181	FLOAT	4-Byte	Humidity to gas sensor Low alarm threshold	R	This is the humidity supplied to gas sensor (CO)
2183	FLOAT	4-Byte	Humidity to gas sensor Hysteresis	R	
2185	INT16U	2-Byte	Relay port/status	R	See explanation at 2136
		Mea	surement Channel Information	on	
2200	INT16U	2-Byte	<b>O2</b> : unit+resolution+type	R	These are the fix information that doesn't need to be read frequently.
2201	INT16U	2-Byte	CO2: unit+resolution+type	R	
2202	INT16U	2-Byte	<b>CO</b> : unit+resolution+type	R	
2203	INT16U	2-Byte	<b>H2O</b> : unit+resolution+type	R	
2204	INT16U	2-Byte	Oil: unit+resolution+type	R	
2205	INT16U	2-Byte	<b>Pressure</b> : unit+resolution+type	R	
2206	INT16U	2-Byte	Temp.: unit+resolution+type	R	
2207	INT16U	2-Byte	Humidity Gas Sensor: unit+resolution+type	R	
		Measur	ement Channel Status and Va	alues	
2300	INT16U	2-Byte	Status (channels 18)	R	Status information,
2301	FLOAT	4-Byte	Channel Value <b>02</b>	R	with 2 bits:

	$\overline{}$	
יכ.	J	Ú

2303	FLOAT	4-Byte	Channel Value CO2	R	
2305	FLOAT	4-Byte	Channel Value <b>CO</b>	R	_
2307	FLOAT	4-Byte	Channel Value <b>H20</b>	R	_
2309	FLOAT	4-Byte	Channel Value <b>Oil</b>	R	_
2311	FLOAT	4-Byte	Channel Value <b>Pressure</b>	R	00 all OK - 01 Out of range
2313	FLOAT	4-Byte	Channel Value <b>Temperature</b>	R	_ 10 Sensor broken
2315	FLOAT	4-Byte	Channel Value <b>Humidity Gas S.</b>	R	11 Other error
		Device sp	ecific System Information for	r <b>S60</b>	5
2600	INT32U	4-Byte	General System status	R	Bit1 bit0(water level): 00: normal 01: low level indicator 10: high level indicator 11: N/A  Bit 2: O2 Bit 3: CO2 Bit 4: CO Bit 5: H2O Bit 6: Oil Bit 7: Pressure Bit 8: temperature Bit 9: Humidity for Gas. 0: normal 1: calibration overdue.
2602	INT32U	4-Byte	Alarm status 1: Bit 0-1: O2 Bit 2-3: CO2 Bit 4-5: CO Bit 6-7: H2O Bit 8-9: Oil Bit 10-11:Pressure Bit 12-13: Temp. Bit 14-15: Humidity Gas Bit 16-32: reserved for future	R	2 bits for alarm status: 00: no alarm 01: High alarm 10: Low alarm 11: N/A In Total 16 alarm status
2604	FLOAT	4-Byte	Total work time	R	Counted from when? The total operating time since delivery
2606	FLOAT	4-Byte	Minimum value O2	R	When are the
2608	FLOAT	4-Byte	Maximum value O2	R	min/max values rest? Need to define
2610	FLOAT	4-Byte	Minimum value CO2	R	

48 \$605



2612	FLOAT	4-Byte	Maximum value CO2	R	this clearly		
2614	FLOAT	4-Byte	Minimum value CO	R			
2616	FLOAT	4-Byte	Maximum value CO	R			
2618	FLOAT	4-Byte	Minimum value H2O	R			
2620	FLOAT	4-Byte	Maximum value H2O	R			
2622	FLOAT	4-Byte	Minimum value Oil	R			
2624	FLOAT	4-Byte	Maximum value Oil	R			
2626	FLOAT	4-Byte	Minimum value pressure	R			
2628	FLOAT	4-Byte	Maximum value pressure	R			
2630	FLOAT	4-Byte	Minimum value temperature	R			
2632	FLOAT	4-Byte	Maximum value temperature	R			
2633	FLOAT	4-Byte	Minimum value humidity gas	R			
2635	FLOAT	4-Byte	Maximum value humidity gas	R			
2637	FLOAT	4-Byte	Minimum value barometer	R			
2639	FLOAT	4-Byte	Maximum value barometer	R			
Specific Settings of integrated Sensors							
10000		BYTE[1 K]	O2 sensor settings	R/W			
10000	INT32U	4-Byte	Serial number	R	Factory S/N		
10002	INT16U	2-Byte	High = 0xFF, Low byte HW version	R	HW version 10 = 1.0		
10003	INT32U	4-Byte	FW version (low 3 Bytes used)	R	New firmware version format: XX.XX.XX		
10005	INT16U	2-Byte	Enable user calibration (0 disable,1 enable)	R/W			
10006	INT16U	2-Byte	User calibration point count	R/W			
10007	FLOAT	4-Byte	User calibration point 0 reference value	R/W			
10009	FLOAT	4-Byte	User calibration point 0 actual value	R/W			
10011	FLOAT	4-Byte	User calibration point 1 reference value (Reserved)	R/W			
10013	FLOAT	4-Byte	User calibration point 1 actual value (Reserved)	R/W			
10500		BYTE[1 K]	CO2 sensor settings	R/W			



10500	INT32U	4-Byte	Serial number	R	Factory S/N
10502	INT16U	2-Byte	High = 0xFF, Low byte HW version	R	HW version 10 = 1.0
10503	INT32U	4-Byte	FW version (low 3 Bytes used)	R	New firmware version format: XX.XX.XX
10505	INT16U	2-Byte	Enable user calibration (0 disable,1 enable)	R/W	
10506	INT16U	2-Byte	User calibration point count	R/W	
10507	FLOAT	4-Byte	User calibration point 0 reference value	R/W	
10509	FLOAT	4-Byte	User calibration point 0 actual value	R/W	
10511	FLOAT	4-Byte	User calibration point 1 reference value (Reserved)	R/W	
10513	FLOAT	4-Byte	User calibration point 1 actual value (Reserved)	R/W	
11000		BYTE[1 K]	CO sensor settings	R/W	
11000	INT32U	4-Byte	Serial number	R	Factory S/N
11002	INT16U	2-Byte	High = 0xFF, Low byte HW version	R	HW version 10 = 1.0
11003	INT32U	4-Byte	FW version (low 3 Bytes used)	R	New firmware version format: XX.XX.XX
11005	INT16U	2-Byte	Enable user calibration (0 disable,1 enable)	R/W	
11006	INT16U	2-Byte	User calibration point count	R/W	
11007	FLOAT	4-Byte	User calibration point 0 reference value	R/W	
11009	FLOAT	4-Byte	User calibration point 0 actual value	R/W	
11011	FLOAT	4-Byte	User calibration point 1 reference value (Reserved)	R/W	
11013	FLOAT	4-Byte	User calibration point 1 actual value (Reserved)	R/W	
11500		BYTE[1 K]	<b>H2O sensor settings</b> (this is the dew point sensor)	R/W	
11500	INT32U	4-Byte	Serial number	R	Factory S/N
11502	INT16U	2-Byte	High = 0xFF, Low byte HW version	R	HW version 10 = 1.0



4-Byte	FW version (low 3 Bytes used)	R	New firmware version format: XX.XX.XX
2-Byte	Alternative dew point unit	R/W	Belongs to H2O
4-Byte	Static pressure (in bar(g))	R/W	sensor
BYTE[1 K]	Oil sensor settings	R/W	
4-Byte	Serial number	R	Factory S/N
2-Byte	High = 0xFF, Low byte HW version	R	HW version 10 = 1.0
4-Byte	FW version (low 3 Bytes used)	R	New firmware version format: XX.XX.XX
BYTE[1 K]	<b>Dew Point settings</b> (this is for gas sensor)	R/W	
4-Byte	Serial number	R	Factory S/N
2-Byte	High = 0xFF, Low byte HW version	R	HW version 10 = 1.0
4-Byte	FW version (low 3 Bytes used)	R	New firmware version format: XX.XX.XX
4-Byte	Atmospheric pressure (in hPa)	R/W	
	2-Byte 4-Byte BYTE[1 K] 4-Byte 2-Byte 4-Byte BYTE[1 K] 4-Byte 2-Byte 4-Byte 4-Byte	2-Byte Alternative dew point unit  4-Byte Static pressure (in bar(g))  BYTE[1 K] Oil sensor settings  4-Byte Serial number  2-Byte High = 0xFF, Low byte HW version  4-Byte FW version (low 3 Bytes used)  BYTE[1 Dew Point settings (this is for gas sensor)  4-Byte Serial number  2-Byte High = 0xFF, Low byte HW version  4-Byte FW version (low 3 Bytes used)  4-Byte FW version (low 3 Bytes used)  4-Byte Atmospheric pressure (in	2-Byte Alternative dew point unit R/W 4-Byte Static pressure (in bar(g)) R/W  BYTE[1 K] Oil sensor settings R/W 4-Byte Serial number R  2-Byte High = 0xFF, Low byte HW version R  4-Byte FW version (low 3 Bytes used) R/W  BYTE[1 Dew Point settings (this is for gas sensor) R/W  4-Byte Serial number R  2-Byte High = 0xFF, Low byte HW version R  4-Byte FW version (low 3 Bytes used) R  4-Byte FW version (low 3 Bytes R  A-Byte FW version (low 3 Bytes R  A-Byte FW version (low 3 Bytes R  A-Byte FW version (low 3 Bytes used) R

### SUTO iTEC GmbH

Grißheimer Weg 21 D-79423 Heitersheim Germany

Tel: +49 (0) 7634 50488 00

Email: <a href="mailto:sales@suto-itec.com">sales@suto-itec.com</a>
Website: <a href="mailto:www.suto-itec.com">www.suto-itec.com</a>

# SUTO iTEC (ASIA) Co., Ltd.

Room 10, 6/F, Block B, Cambridge Plaza 188 San Wan Road, Sheung Shui, N.T. Hong Kong

Tel: +852 2328 9782

Email: <a href="mailto:sales.asia@suto-itec.com">sales.asia@suto-itec.com</a>
Website: <a href="mailto:www.suto-itec.com">www.suto-itec.com</a>

## SUTO iTEC Inc.

5460 33rd St SE Grand Rapids, MI 49512 USA

Tel: +1 (616) 800-7886

Email: <a href="mailto:sales.us@suto-itec.com">sales.us@suto-itec.com</a>
Website: <a href="mailto:www.suto-itec.com">www.suto-itec.com</a>

All rights reserved ©

Modifications and errors reserved

S605\_im\_en\_2023-2.odt